

**LOOSEBOX / LYING AREA**

**BELMONDO® Trend**

**For extra lying comfort – horse-friendly softness**

- heat insulating air cushion for horse-friendly softness
- sealing lips on lower side for better hygiene
- for safe getting up and lying down
- appropriate step absorption relieves joints
- less noise, less dust and ammonia exposure contribute to a healthier climate in the stable
- material impermeable to water and resistant to urine
- minimizes the amount of bedding, costs and labour



|                             |  |
|-----------------------------|--|
| <b>Surface:</b>             | horseshoe profile with wear resistant top coat |
| <b>Lower side:</b>          | air cushions with sealing lips                 |
| <b>Material:</b>            | rubber   |
| <b>Thickness:</b>           | 28 mm  |
| <b>Dimensions:</b>          | 1 m x 1 m                                      |
| <b>Weight:</b>              | approx. 22 kg                                  |
| <b>Puzzle shaped edges:</b> | on 4 sides                                     |



**horseshoe profile with top coat**



**air cushions with sealing lips on lower side**

Special dimensions are generally available on request. All BELMONDO® products are suitable for looseboxes as well as open stables.

To bind fluids (urine) a small amount of absorbent bedding is always required in the loosebox. We also recommend a minimal amount of bedding in the open stable.



**Laying:** according to BELMONDO® laying instructions

**Tolerances:** thickness: +3/-2 mm  
length/width: +/- 1.5 % (DIN ISO 3302-1 tolerance class M4)  
The tolerances also apply within the mat.

The products presented here should only be used for the listed applications. Technical specifications are subject to alterations. All sales are subject to our General Terms and Conditions. Legal basis is the German version of the document.